

MOGAT2 Antibody



PACO10544

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

MOGAT2

Uniprot

Q3SYC2

Synonyms:

monoacylglycerol O-acyltransferase 2; MOGAT2; DGAT2L5; FLJ22644; MGAT2; MGC119183; MGC119184; MGC119185 ;

Immunogen:

Human MOGAT2.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A