

MORC1 Antibody



PACO10546

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

MORC1

Uniprot

Q86VD1

Synonyms:

MORC family CW-type zinc finger 1; MORC1; CT33; MORC; ZCW6 ;

Immunogen:

Human MORC1.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A