

MORF4L1 Antibody



PACO10547

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

MORF4L1

Uniprot

Q9UBU8

Synonyms:

mortality factor 4 like 1; MORF4L1; Eaf3; FWP006; HsT17725; MGC10631; MORFRG15; MRG15; S863-6 ;

Immunogen:

Human MORF4L1.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A