

MPP5 Antibody



PACO10556

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

MPP5

Uniprot

Q8N3R9

Synonyms:

membrane protein, palmitoylated 5 (MAGUK p55 subfamily member 5); MPP5; FLJ12615; PALS1 ;

Immunogen:

Human MPP5.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A