

# MPPED1 Antibody



PACO10559

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC, IF

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

MPPED1

**Uniprot**

O15442

**Synonyms:**

metallophosphoesterase domain containing 1; MPPED1;239AB; C22orf1; FAM1A; FLJ50009; FLJ54633; FLJ59965; FLJ78907; MGC88045 ;

**Immunogen:**

Human MPPED1.

**Storage:**

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A