

MPPED2 Antibody



PACO10560

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:

Protein Background:

Gene ID:

MPPED2

Uniprot

Q15777

Synonyms:

metallophosphoesterase domain containing 2; MPPED2;239FB; C11orf8 ;

Immunogen:

Human MPPED2.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A