

MPRIP Antibody



PACO10561

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

Protein Background:

Gene ID:

MPRIP

Uniprot

Q6WCQ1

Synonyms:

myosin phosphatase Rho interacting protein; MPRIP; KIAA0864; M-RIP; RHOIP3; p116Rip ;

Immunogen:

Human MPRIP.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A