

# MRAP Antibody



PACO10565

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

MRAP

**Uniprot**

Q8TCY5

**Synonyms:**

melanocortin 2 receptor accessory protein; MRAP; B27; C21orf61; FALP; FGD2; GCCD2 ;

**Immunogen:**

Human MRAP.

**Storage:**

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A