

MRPL19 Antibody



PACO10577

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

MRPL19

Uniprot

P49406

Synonyms:

mitochondrial ribosomal protein L19; MRPL19; KIAA0104; L19mt; MGC20675; MRP-L15; MRP-L19; MRPL15; RLX1; RPML15 ;

Immunogen:

Human MRPL19.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A