

MRPL24 Antibody



PACO10581

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:**Protein Background:****Gene ID:**

MRPL24

Uniprot

Q96A35

Synonyms:

mitochondrial ribosomal protein L24; MRPL24; FLJ20917; L24mt; MGC22737; MGC9831; MRP-L18; MRP-L24 ;

Immunogen:

Human MRPL24.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A