

MRPL36 Antibody



PACO10586

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

MRPL36

Uniprot

Q9P0J6

Synonyms:

mitochondrial ribosomal protein L36; MRPL36; BRIP1; L36mt; MGC104245; MRP-L36; PRPL36; RPMJ ;

Immunogen:

Human MRPL36.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A