

MRPL49 Antibody



PACO10593

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

MRPL49

Uniprot

Q13405

Synonyms:

mitochondrial ribosomal protein L49; MRPL49; C11orf4; L49mt; MGC10656; NOF; NOF1 ;

Immunogen:

Human MRPL49.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A