

MRPS10 Antibody



PACO10596

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

MRPS10

Uniprot

P82664

Synonyms:

mitochondrial ribosomal protein S10; MRPS10; FLJ10567; MRP-S10; PNAS-122 ;

Immunogen:

Human MRPS10.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A