

MRPS11 Antibody



PACO10597

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

MRPS11

Uniprot

P82912

Synonyms:

mitochondrial ribosomal protein S11; MRPS11; FLJ22512; FLJ23406; HCC-2 ;

Immunogen:

Human MRPS11.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A