

# MRPS21 Antibody



PACO10605

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IHC

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

MRPS21

**Uniprot**

P82921

**Synonyms:**

mitochondrial ribosomal protein S21; MRPS21; MDS016; MRP-S21; RPMS21 ;

**Immunogen:**

Human MRPS21.

**Storage:**

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A