

MRPS23 Antibody



PACO10607

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:

Protein Background:

Gene ID:

MRPS23

Uniprot

Q9Y3D9

Synonyms:

mitochondrial ribosomal protein S23; MRPS23; CGI-138; HSPC329; MRP-S23 ;

Immunogen:

Human MRPS23.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A