

MRPS28 Antibody



PACO10612

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

MRPS28

Uniprot

Q9Y2Q9

Synonyms:

mitochondrial ribosomal protein S28; MRPS28; FLJ22853; HSPC007; MRP-S28; MRP-S35; MRPS35 ;

Immunogen:

Human MRPS28.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A