

MRPS35 Antibody



PACO10615

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:

Protein Background:

Gene ID:

MRPS35

Uniprot

P82673

Synonyms:

mitochondrial ribosomal protein S35; MRPS35; DKFZp762P093; HDCMD11P; MDS023; MGC104278; MRP-S28; MRPS28 ;

Immunogen:

Human MRPS35.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A