

MRTO4 Antibody



PACO10618

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

MRTO4

Uniprot

Q9UKD2

Synonyms:

mRNA turnover 4 homolog (S. cerevisiae); MRTO4; C1orf33; MRT4; dJ657E11.4 ;

Immunogen:

Human MRTO4.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A