

MST1R Antibody



PACO10629

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

MST1R

Uniprot

Q04912

Synonyms:

macrophage stimulating 1 receptor (c-met-related tyrosine kinase); MST1R; CD136; CDw136; PTK8; RON ;

Immunogen:

Human MST1R.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A