

MTERF1 Antibody



PACO10635

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

MTERF1

Uniprot

Q99551

Synonyms:

mitochondrial transcription termination factor; MTERF; MGC131634 ;

Immunogen:

Human MTERF.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A