

# MTG1 Antibody



PACO10638

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

MTG1

**Uniprot**

Q9BT17

**Synonyms:**

mitochondrial GTPase 1 homolog (S. cerevisiae); MTG1; GTP; GTPBP7; RP11-108K14.2 ;

**Immunogen:**

Human MTG1.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A