

MTHFS Antibody



PACO10641

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

MTHFS

Uniprot

P49914

Synonyms:

5,10-methenyltetrahydrofolate synthetase (5-formyltetrahydrofolate cyclo-ligase);
MTHFS; FLJ30410; HsT19268 ;

Immunogen:

Human MTHFS.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A