

MTMR8 Antibody



PACO10650

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

MTMR8

Uniprot

Q96EF0

Synonyms:

myotubularin related protein 8; MTMR8; FLJ20126 ;

Immunogen:

Human MTMR8.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A