

MTRF1 Antibody



PACO10653

Product Information

Size:

50ul

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

MTRF1

Uniprot

O75570

Synonyms:

mitochondrial translational release factor 1; MTRF1; MGC47721; MRF1; MTTRF1; RF1 ;

Immunogen:

Human MTRF1.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A