

MTUS1 Antibody



PACO10654

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

MTUS1

Uniprot

Q9ULD2

Synonyms:

mitochondrial tumor suppressor 1; MTUS1; ATIP; DKFZp586D1519; DKFZp686F20243; FLJ14295; KIAA1288; MP44; MTSG1 ;

Immunogen:

Human MTUS1.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A