

# MYBPC3 Antibody



PACO10671

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:****Protein Background:****Gene ID:**

MYBPC3

**Uniprot**

Q14896

**Synonyms:**

myosin binding protein C, cardiac; MYBPC3; CMH4; DKFZp779E1762; FHC; MYBP-C ;

**Immunogen:**

Human MYBPC3.

**Storage:**

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A