

MYCBP Antibody



PACO10673

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

Protein Background:

Gene ID:

MYCBP

Uniprot

Q99417

Synonyms:

c-myc binding protein; MYCBP; AMY-1 ;

Immunogen:

Human MYCBP.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A