

MYCBPAP Antibody



PACO10674

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

MYCBPAP

Uniprot

Q8TBZ2

Synonyms:

MYCBP associated protein; MYCBPAP; AMAP-1; AMAP1; DKFZp434N1415 ;

Immunogen:

Human MYCBPAP.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A