

PACO10677

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

MYEOV

Uniprot

Q96EZ4

Synonyms:

myeloma overexpressed (in a subset of t(11;14) positive multiple myelomas); MYEOV; OCIM ;

Immunogen:

Human MYEOV.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A