

MYH2 Antibody



PACO10682

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

MYH2

Uniprot

Q9UKX2

Synonyms:

myosin, heavy chain 2, skeletal muscle, adult; MYH2; MYH2A; MYHSA2; MYHas8; MyHC-2A; MyHC-IIa ;

Immunogen:

Human MYH2.

Storage:

0.1M NaHCO₃, 0.1M Glycine, 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C; C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A