

# MYH9 Antibody



PACO10683

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC, IF

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

MYH9

**Uniprot**

P35579

**Synonyms:**

myosin, heavy chain 9, non-muscle; MYH9; DFNA17; EPSTS; FTNS; MGC104539; MHA; NMHC-II-A; NMMHCA ;

**Immunogen:**

Human MYH9.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A