

MYL3 Antibody



PACO10687

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

MYL3

Uniprot

P08590

Synonyms:

myosin, light chain 3, alkali; ventricular, skeletal, slow; MYL3; CMH8; MLC1SB; MLC1V; VLC1 ;

Immunogen:

Human MYL3.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A