

MYL7 Antibody



PACO10690

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

MYL7

Uniprot

Q01449

Synonyms:

myosin, light chain 7, regulatory; MYL7; MYL2A; MYLC2A ;

Immunogen:

Human MYL7.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A