

MYLK2 Antibody



PACO10694

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

MYLK2

Uniprot

Q9H1R3

Synonyms:

myosin light chain kinase 2; MYLK2; KMLC; MLCK; MLCK2; skMLCK ;

Immunogen:

Human MYLK2.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A