

MYLK3 Antibody



PACO10695

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

MYLK3

Uniprot

Q32MK0

Synonyms:

myosin light chain kinase 3; MYLK3; MGC126319; MGC126320; MLCK; MLCK2; caMLCK ;

Immunogen:

Human MYLK3.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A