

# MYLPF Antibody



PACO10696

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

MYLPF

**Uniprot**

Q96A32

**Synonyms:**

myosin light chain, phosphorylatable, fast skeletal muscle; MYLPF; DKFZp779C0757; HUMMLC2B; MGC13450; MRLC2; MYL11 ;

**Immunogen:**

Human MYLPF.

**Storage:**

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A