

# MYO16 Antibody



PACO10697

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

MYO16

**Uniprot**

Q9Y6X6

**Synonyms:**

myosin XVI; MYO16; KIAA0865; MYR8; Myo16b ;

**Immunogen:**

Human MYO16.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A