

MYO1A Antibody



PACO10699

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

Protein Background:

Gene ID:

MYO1A

Uniprot

Q9UBC5

Synonyms:

myosin IA; MYO1A; BBMI; DFNA48; MIHC; MYHL ;

Immunogen:

Human MYO1A.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A