

MYO1D Antibody



PACO10701

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:**Protein Background:****Gene ID:**

MYO1D

Uniprot

O94832

Synonyms:

myosin ID; MYO1D; KIAA0727; myr4 ;

Immunogen:

Human MYO1D.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A