

MYO1E Antibody



PACO10702

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

MYO1E

Uniprot

Q12965

Synonyms:

myosin IE; MYO1E; HuncM-IC; MGC104638; MYO1C ;

Immunogen:

Human MYO1E.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A