

# MYO7A Antibody



PACO10703

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IF

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

MYO7A

**Uniprot**

Q13402

**Synonyms:**

myosin VIIA; MYO7A; DFNA11; DFNB2; MYOVIIA; MYU7A; NSRD2; USH1B ;

**Immunogen:**

Human MYO7A.

**Storage:**

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A