

# MYO7B Antibody



PACO10704

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IHC

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

MYO7B

**Uniprot**

Q6PIF6

**Synonyms:**

myosin VIIB; MYO7B; DKFZp686A08248 ;

**Immunogen:**

Human MYO7B.

**Storage:**

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A