

MYOC Antibody



PACO10706

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

MYOC

Uniprot

Q99972

Synonyms:

myocilin, trabecular meshwork inducible glucocorticoid response; MYOC; GLC1A; GPOA; JOAG; JOAG1; TIGR; myocilin ;

Immunogen:

Human MYOC.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A