

MYOCD Antibody



PACO10707

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

Protein Background:

Gene ID:

MYOCD

Uniprot

Q8IZQ8

Synonyms:

myocardin; Myocd; BSAC2A; Srfcp ;

Immunogen:

Human Myocd.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A