

# N6AMT2 Antibody



PACO10715

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:****Protein Background:****Gene ID:**

EEF1AKMT1

**Uniprot**

Q8WVE0

**Synonyms:**

N-6 adenine-specific DNA methyltransferase 2 (putative); N6AMT2; ESP13 ;

**Immunogen:**

Human N6AMT2.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A