

# NADK Antibody



PACO10717

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC, IF

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

NADK

**Uniprot**

O95544

**Synonyms:**

NAD kinase; NADK; FLJ13052; FLJ37724; dJ283E3.1 ;

**Immunogen:**

Human NADK.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A