

NARF Antibody



PACO10732

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

NARF

Uniprot

Q9UHQ1

Synonyms:

nuclear prelamin A recognition factor; NARF; DKFZp434G0420; FLJ10067; IOP2 ;

Immunogen:

Human NARF.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A