

NARS2 Antibody



PACO10735

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

NARS2

Uniprot

Q96I59

Synonyms:

asparaginyl-tRNA synthetase 2, mitochondrial (putative); NARS2; FLJ23441; SLM5 ;

Immunogen:

Human NARS2.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A