

# NAT1 Antibody



PACO10737

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

NAT1

**Uniprot**

P18440

**Synonyms:**

N-acetyltransferase 1 (arylamine N-acetyltransferase); NAT1; AAC1; NATI ;

**Immunogen:**

Human NAT1.

**Storage:**

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A